

## Iris 15<sup>™</sup> Scientific CMOS Camera Datasheet

Scientific CMOS, EMCCD and CCD Cameras



## Large Field of View Scientific CMOS

Iris 15 Scientific CMOS camera is designed with a large field of view for live cell microscopy applications. Iris 15 delivers a 15 Megapixel sensor with a 25 millimeter field of view ensuring it can maximize the number of cells captured in a single frame. Iris 15 is equipped with an imaging area that is 50% larger than current sCMOS devices.

The 4.25µm pixels provide highly detailed images across the imaging plane and allow for Nyquist spatial sampling at 40X magnification. The camera can capture dynamic cellular events at 30 frames per second for the full frame, and at thousands of frames per second with regions.

The Iris 15 has a high quantum efficiency and low noise levels to maximize dim signal detection and allowing for the use of shorter exposure times to minimize cellular photo-damage.

The Iris 15 is the ideal camera to leverage the larger fields of view of newer microscopes and to image with high frame rates with high resolution for live-cell microscopy applications.



Primary applications: Light Sheet Microscopy Multi-Color Fluorescence High Content Screening Tiling Applications

- 15 Megapixel (5056 x 2960)
- > 25mm Field of View
- 4.25µm x 4.25µm Pixel Area
- > 30 Frames per Second
- >73% Quantum Efficiency

Features	Advantages
15 Megapixel Sensor 25mm Field of View	Maximize the imaging area and increase the number of cells acquired per frame
Optimized 4.25µm Pixel Size	Image the finer details of your samples and maintain proper spatial sampling at 40X magnification
Fast Frame Rates	Capture dynamic cellular events with high temporal resolution
High Quantum Efficiency	Maximizes the ability to detect signals, enables shorter exposure times and minimizes photo-toxicity
Low Read Noise	Maximize ability to detect extremely faint fluorescence signals
Large Dynamic Range	Measure both bright and dim signal levels within the same image
Programmable Scan Mode	Easily synchronize and control acquisitions with the rolling shutter readout



Specifications	Camera Performance
Sensor	Gpixel GSense 5130 Scientific CMOS sensor
Active Array Size	5056 x 2960 (15 Megapixel)
Pixel Area	4.25μm x 4.25μm (18.06μm2)
Sensor Area	21.49mm x 12.61mm 24.9mm diagonal
Peak QE%	>73%
Read Noise	1.5e-
Full-Well Capacity	13,000e-
Bit Depth	16-bit
Readout Mode	Rolling Shutter Effective Global Shutter Programmable San Mode (PCI-E only)
Binning	2x2 (on FPGA)


Cooling Performance	Sensor Temperature	Dark Current
Air Cooled	0°C @ 30°C Ambient	0.5e-/pixel/second

Specification	Camera Interface	
Digital Interface	PCle, USB 3.0	
Lens Interface	F-Mount	
Mounting Points	1/4"-20 mounting point on each side	

Triggering Mode	Function
Input Trigger Medes	Trigger-First: Sequence triggered on first rising edge
Input Trigger Modes	Edge: Each frame triggered on rising edge
Output Trigger Modes	Any Row: Expose signal is high while any row is acquiring data
	All Rows: Effective Global Shutter – Expose signal is high when all rows are acquiring data
	Signal is high for set Exposure time
	Rolling Shutter: Effective Global Shutter – Expose signal is high when all rows are acquiring data
	Signal is High for set Exposure time – Readout Time
	Line Output: Expose signal provides rising edge for each row advanced by the rolling shutter
	readout
Output Trigger Signals	Expose Out, Read Out, Trigger Ready


Programmable Scan Mode	Function
	Auto: Normal camera operation
Scan Modes	Line Delay: Control rolling shutter propagation rate by adding delays to the line time
	Scan Width: Control number of rows between reset and readout signal in the rolling shutter
	Down: Rolling shutter readout begins at the top of the sensor
Coop Direction	Up: Rolling shutter readout begins at the bottom of the sensor
Scan Direction	Down/Up Alternate: Rolling shutter readout alternates direction after starting at the top of the
	sensor








Frame Rate			
PCI-Express	USB 3.0		
16-bit	16-bit		
30	10		
59	31		
174	91		
695	321		
	PCI-Express   16-bit   30   59   174		

Accessories (Included)		
PCIe Card/Cable	Manual	
Trigger Cable	QuickStart Guide	
Power Supply		









Teledyne Photometrics is a registered trademark. Iris 15 is a trademark of Teledyne Photometrics. All other brand and product names are the trademarks of their respective owners.

## Specifications in this datasheet are subject to change.

Refer to the Teledyne Photometrics website for most current specifications.

www.photometrics.com info@photometrics.com / tel: +1 520.889.9933

