

QIMAGING QUICK INSTALLATION GUIDE

Install the Camera-Interface Software

1. Install the IEEE-1394 card in your computer (refer to the video <http://www.qimaging.com/qitube/?v=installing-pci-card>). **Note:** Please use the power supply which came with your camera.
2. Power on the computer and insert the QImaging USB stick into your USB drive.
3. Install camera drivers by choosing "QI_setup.exe".
4. Follow the installation instructions on screen.
5. Restart your computer and connect the camera to the PC.

EXi™ Blue





Retiga™ 2000DC & 4000DC



Rolera™ EM-C²



1. Plug one end of the cable into one of the camera's Firewire sockets (either socket is fine). Plug the other end of the cable into a Firewire port of your computer.



2. Connect the power cord to the power supply, plug the power cord into a wall outlet & connect the power supply to the camera.



3. Power on your computer. Then power on the camera by pushing the power button.

Note: The Retiga 4000DC camera pictured above is displayed with a standard c-mount adapter. The camera, however, ships with an f-mount adapter attached and a c-mount adapter as an accessory. These adapters are interchangeable and either one can be used depending on the adapter requirements of the microscope.

More detailed installation instruction can be found in the camera instruction manual on the USB stick.

Contact QImaging Customer Service for additional information: 1-800-874-9789 or www.qimaging.com



www.qimaging.com

©2015 QImaging. All rights reserved.

03-0025 Rev A2