



**TELEDYNE
PHOTOMETRICS**
Everywhereyoulook™



Prime 95B™ Scientific CMOS Camera

CUSTOMER REFERENCE

Spectroscopy Application In High Flux Plasma Generator Masnum-PSI

Dr. Hennie van der Meiden and Mr. Gijs Akkermans

Dutch Institute for Fundamental Energy Research (DIFFER), Eindhoven, The Netherlands.

BACKGROUND

With the transition to sustainable energy and new energy technologies, research into controlled nuclear fusion is being carried out. In a worldwide collaboration, the tokamak ITER is being developed and built in Cadarache, France. This is one of the first steps for realizing power plants that can generate clean and reliable energy for the future.

The principle design concept of ITER is a doughnut-shaped vessel where the fusion fuel, a hot plasma of hydrogen isotopes, is contained by high magnetic fields and fusion power can be produced at a plasma temperature of ~100 million °C. Fusion products like helium and impurities have to be exhausted from the core plasma which results in high power loads of 10 MW/m² (continuous), peaking at over 1 GW/m² on plasma facing components (PFCs) inside the exhaust systems. To address this challenge, a linear plasma generator Magnum-PSI was built to study the plasma-wall interactions (Fig.1).

(continued...)




Figure 1: A side cutaway view of the linear plasma generator Magnum-PSI. The plasma (magenta-colored beam) exhausted by the plasma source (right side) is confined by the max 2.5 Tesla magnetic field (super conductive magnet (dark blue)) and deposited on a sample (power flux <math><40 \text{ MW/m}^2</math> continuously and transients >1GW/m²). Different diagnostics are used for measuring the plasma conditions.

BACKGROUND

PhD student Gijs Akkermans and Dr Hennie van der Meiden are involved in the research with the Magnum-PSI in order to test different materials under different plasma conditions similar to those in ITER and beyond. As he told us, “the power load in some areas of the tokamak is too high and can damage the plasma-facing components”.

In their experiments, a continuous or pulsed beam of plasma is focused onto a sample. Different diagnostics are used to monitor the plasma conditions such as temperature, density, and calorimetry. Among these, a spectrometer system is used for measuring the population of states of atoms and molecules by collecting the photons originating from these particles. Inside the spectrometer, these photons are separated by wavelength and imaged using a scientific camera.

“With the Prime 95B 25mm I can get a 10 μ s exposure time with the fast mode, letting us image the induced plasma pulses.”

CHALLENGE

The cameras used for the spectroscopy need to be highly UV sensitive as the wavelength of interest is 340-500 nm, with good signal-to-noise ratio for UV wavelengths.

For experiments with continuous plasma, different exposure times are needed depending on the intensity and wavelength of the signal, making the dynamic range important. The plasma pulses have a duration of 0.5 ms, and need to be sampled with time scales shorter than in 50 μ s, making very fast shutter times desirable. In addition both types of experiment need to capture the entire field of view without significant magnification, meaning that cameras with a large field of view are necessary.

SOLUTION

The Prime 95B 25mm delivers on each of the challenges presented by spectroscopy and imaging of the radiated plasma photons, offering great quantum efficiency and low noise at UV wavelengths. As mentioned by Dr. Hennie van der Meiden “with the Prime 95B 25mm I can get a 10 μ s exposure time with the fast mode, letting us image the induced plasma pulses and plasma fluctuations.”

The large sensor and wide field of view offered by the Prime 95B 25mm ensures that a large area can be imaged without having to compromise high speed and high resolution, all with very low levels of readout noise and dark current.





Figure 2: The left figure shows a spectrum at middle fiber low resolution setting. The right figure shows light from the plasma beam imaged on a fiber array and the spectra imaged with the Prime 95B. Each fiber corresponds to a certain plasma position of the beam, i.e. a radial profile of the beam is measured. The highest peak corresponds to a Balmer beta line (488 nm, $n=4 > n=2$).

