



Kinetix Scientific CMOS Camera

**CUSTOMER REFERENCE**

## Functional Calcium Imaging Using Light Sheet Microscopy

**Dr Asaph Zylbertal and Dr Isaac Bianco**

Department of Neuroscience, Physiology &amp; Pharmacology, University College London, UK

**BACKGROUND**

The Bianco Lab at UCL aims to understand the logic of brain circuits that control behaviour. Using two photon and light sheet microscopy with small, optically transparent larval zebrafish expressing genetically encoded calcium indicators (GECIs), researchers such as Dr. Asaph Zylbertal record activity at single-cell resolution throughout the brain.

Dr. Zylbertal explained more about his work, "I am using our light sheet microscope to look at tens of thousands of neurons expressing GECIs across the whole zebrafish brain, to see activity and estimate what each individual neuron in the fish brain is doing at any given time."

"After acquiring the raw data, we use a segmentation algorithm to find neurons in each imaging plane, we then deconvolve the calcium signal from these neurons to estimate spikes, this gives us a representation of activity across the brain."

While some light sheet experiments look at the anatomy of the zebrafish, the aim of Dr. Zylbertal is to observe functional activity within the zebrafish brain.

**“The Kinetix allows us to imaging a large field of view at high speeds and with high sensitivity which is ideal for our research.”**

**CHALLENGE**

Calcium imaging requires high speeds in order to capture the dynamic signalling between neurons. However, to monitor populations of cells throughout the brain, it is also necessary to capture large imaging volumes. Thus, the detector needs to be capable of high sensitivity imaging and both high speed and a large field of view (FOV) for this high throughput application.

Dr. Zylbertal told us about the experimental challenges, "When we are using fast indicators, we would like to be able to image as many volumes as possible over time, all with high resolution so we can integrate lots of pixels in order to get a higher signal to noise."

SOLUTION

The Kinetix is an ideal solution for this application, featuring extremely high speeds (500 fps across full frame), high sensitivity (<1 e- read noise and 95% QE signal collection) and a huge FOV (29 mm diagonal), all in one camera.


The Kinetix Speed Mode operates at 500 Hz across the full 10 megapixel sensor, enabling the Kinetix to capture rapid calcium signals from neurons in real time. Combining this with the high sensitivity and large FOV allows for volumetric functional imaging across the whole zebrafish brain, while retaining high resolution to resolve individual neurons.

Dr. Zylbertal outlined his experience with the Kinetix, "The advanced triggering is very important. We synchronise the camera with the light sheet and use pseudo-global shutter"

"With the Kinetix we are imaging at a larger FOV and a higher speed than previously, which is what we need for this application. Experiments can last for hours and we want to react to the images as they come. With the Kinetix we have the higher signal to noise ratio needed to extract activity from the neurons."

REFERENCE

Wolf, S., Dubreuil, A.M., Bertoni, T. et al. (2017) *Sensorimotor computation underlying phototaxis in zebrafish*. Nature Communications 8, 651. <https://doi.org/10.1038/s41467-017-00310-3>



**Figure 1:** Image of the brain of a zebrafish larva, taken with the Kinetix sCMOS. This sample is expressing the GCaMP6f calcium indicator under the elavl3 promotor (elavl3:GCaMP6f, Wolf *et al.* 2017). Dorsal view, view is approx. 420 x 710 x 150  $\mu\text{m}^3$