





Large Field of View BSI Scientific CMOS

Discovery depends on every photon

The Prime 95B 25 mm delivers an extremely large imaging area and combines this with the near-perfect sensitivity of a Backside Illuminated (BSI) Scientific CMOS sensor. The Prime 95B 25 mm sensor converts up to 95% of incident photons into measurable signal, and provides an unmatched 25 mm field of view. The Prime 95B 25 mm is optimally positioned to maximize detection and imaging throughput when mounted on the Nikon ECLIPSE Ti2 microscope.

The extreme sensitivity not only allows fainter signals to be detected, it provides the flexibility to increase frame rates, or turn down the excitation intensity to reduce cellular photo-damage. The Prime 95B 25 mm improves the field of view and maintains the high frame rates and extremely low read noise that has made sCMOS so popular for live-cell imaging.

- ▶ 95% Quantum Efficiency
- > 25mm Field of View (1608x1608)
- ▶ 1.6e- Read Noise (median)
- > 30fps @ 16-bit / 60 fps @ 12-bit





| Features | Advantages |
|--|---|
| High Quantum Efficiency 95% Peak QE | Maximizes ability to detect weak signals, enables short exposure times for high frame rates, minimizes phototoxicity across a wide range of wavelengths |
| Large 25 mm Field of View | Maximize imaging area and increase throughput |
| Large 11 µm Pixel Size | Maximize light collection while maintaining proper spatial sampling |
| Extremely Low Read Noise | Maximize your ability to detect faint fluorescence |
| Fast Frame Rates | Capture highly dynamic events with high temporal resolution |
| Enhanced Dynamic Range | Measure both bright and dim signal levels within the same image 50,000:1 Dynamic Range (94 dB) |
| Multiple Expose Out Triggering | Control up to four light sources for multi-wavelength acquisitions |
| SMART Streaming | Faster acquisition rates with variable exposures, ideal for multi-probed live cell imaging Compatible with Multiple Expose Out Triggering |





2.6 Megapixel BSI CMOS Sensor

Backside Illuminated Sensor 1.6 e- Read Noise (Median) >95% peak QE 80,000 e- full well 11 x 11 µm pixels 25 mm diagonal

Easily Mounted and Secured

F-Mount Two $\frac{1}{4}$ "-20 mounting holes per side

Convenient Interfaces

16-bit Data

• 30 fps

12-bit Data

• 60 fps

Multiple Cooling Options

Forced Air Cooling

- -20°C Cooling
- Selectable Fan Speed

Liquid Cooling

- -25°C Cooling
- Leak-proof, quick-disconnect ports



Advanced Application Triggers

Effective Global Shutter
Up to four selectable expose-out lines




| Specifications | Camera Performance |
|--------------------|--|
| Sensor | GPixel GSense 400 BSI CMOS Gen IV, Grade 1 in imaging area |
| Active Array Size | 1608 x 1608 (2.58 Megapixel) |
| Pixel Area | 11 μm x 11 μm (121 μm²) |
| Sensor Area | 17.69 mm x 17.69 mm 25 mm diagonal |
| Peak QE% | >95% |
| Read Noise | 1.6 e ⁻ (Median) 1.8 e ⁻ (RMS) |
| Full-Well Capacity | 80,000 e- |
| Dynamic Range | 50,000:1 |
| Bit Depth | 16-bit, 12-bit |
| Readout Mode | Rolling Shutter Effective Global Shutter |
| Binning | 2x2 (on FPGA) |
| Linearity | >99.5% |

| Cooling Performance | Sensor Temperature Dark Current | | |
|---------------------|---|--|--|
| Air Cooled | -20°C @ 25°C Ambient | PC Ambient 0.55 e ⁻ /pixel/second | |
| Liquid Cooled | -25°C @ 25°C Ambient 0.3 e ⁻ /pixel/second | | |

| Specification | Camera Performance | |
|-------------------|---|--|
| Digital Interface | PCle | |
| Lens Interface | F-Mount | |
| Mounting Points | 2x 1/4"-20 mounting points per side to prevent rotation | |
| Liquid Cooling | Quick Disconnect Ports | |

| Triggering Mode | Function |
|------------------------|---|
| | Trigger First: Sequence triggered on first rising edge |
| Input Trigger Modes | Edge: Each frame triggered on rising edge |
| | SMART Streaming: Fast iteration through multiple exposure times |
| | First Row: Expose signal is high while first row is acquiring data |
| Output Trigger Modes | Any Row: Expose signal is high while any row is acquiring data |
| | All Rows: Effective Global Shutter – Expose signal is high when all rows are acquiring data |
| Output Trigger Signals | Expose Out (up to four signals), Read Out, Trigger Ready |





| Frame Rate (PCIe interface) | | | | | |
|-----------------------------|--------|--------|--|--|--|
| Array Size | 16-bit | 12-bit | | | |
| 1608 x 1608 | 30 | 60 | | | |
| 1608 x 1200 | 40 | 80 | | | |
| 1608 x 512 | 94 | 188 | | | |
| 1608 x 256 | 188 | 374 | | | |
| 1608 x 128 | 374 | 737 | | | |

Accessories (Included)


PCle Card/Cable Power Supply


USB 3.0 Cable Manuals and QuickStart Guide

Trigger Cables Performance and Gain Calibration Test Data

Accessories (Additional)

Liquid Circulator
Liquid Cooling Tubes







Teledyne Photometrics is a registered trademark.

Prime 95B 25mm is a trademark of Teledyne Photometrics.

All other brand and product names are the trademarks of their respective owners.

Specifications in this datasheet are subject to change. Refer to the Teledyne Photometrics website for most current specifications.



