




1392 x 1040 imaging array
4.65 x 4.65- μm pixels

The CoolSNAP™ cf2 Monochrome camera from Photometrics® incorporates low-noise electronics and moderate CCD cooling to achieve good low-light sensitivity. A megapixel sensor with small, square elements ensures that each image shows extraordinary detail. This feature, along with a high-speed digitizer, shutterless operation, and an interline-transfer CCD, makes the CoolSNAP cf2 Monochrome camera ideal for high-resolution life science imaging applications.



Primary applications
Fixed-cell fluorescence
Pathology
Histology
DIC phase-contrast imaging
Darkfield imaging

Features	Benefits
20-MHz readout	Fast image readout for high-speed focus and image capture
1392 x 1040 imaging array 4.65 x 4.65- μm pixels	Resolves fine detail
Interline-transfer, progressive-scan CCD	Full resolution in every frame
Flexible binning and readout	Increases signal-to-noise performance while increasing the frame rate
IEEE-1394a or PCI interface	High-bandwidth, uninterrupted data transfer with no dropped frames
12-bit digitization	Quantifies bright and dim signals in the same image
Thermoelectric cooling	Low dark current allows longer integration times
C-mount	Easily attaches to microscopes, standard lenses, or optical equipment
Subcompact, fanless design	Low profile allows easy integration
Acquisition software	Captures, analyzes, and saves high-resolution images
PVCam®	Supported by numerous third-party software packages
Circular buffers	Real-time focus
Device sequencing	Precise integration with shutters, filter wheels, etc.
IEEE-1394a compatibility	Windows® XP/Vista 32 and Mac OS X
PCI compatibility	Windows XP/Vista 32, Mac OS X, and Linux® (kernel versions 2.4 and 2.6.8)



Binning	Region		
	1392 x 1040	512 x 512	256 x 256
1 x 1	10	20	38
2 x 2	20	37	62
3 x 3	28	49	77
4 x 4	35	59	90
8 x 8	55	86	114

(Frames per second)

Note: Frame rates are measured at 20 MHz with 0-millisecond exposure times.

Specifications	
CCD image sensor	Sony® ICX205AL; interline-transfer, progressive-scan device with microlenses
CCD format	1392 x 1040 imaging array 4.65 x 4.65-µm pixels 6.5 x 4.8-mm imaging area (optically centered) 1/2" format
Grade	Sony Grade 0
System gain	3 e-/ADU
Linear full well	10,200 e-
Read noise	10 e- rms @ 20 MHz
Nonlinearity	<4%
Digitizer type	12 bits @ 20 MHz
Frame readout	96 ms/frame
CCD temperature	5°C below ambient
Dark current	<1 e-/p/s
Operating environment	15 to 30°C ambient, 0 to 80% relative humidity noncondensing
Dimensions	4.5" x 5.0" x 2.5" (1.9 lbs)
I/O	TTL output while exposing (BNC connector)

Note: Specifications are typical and subject to change.

CoolSNAP is a trademark of Photometrics. Photometrics and PVCam are registered trademarks of Photometrics. Mac OS is a trademark of Apple Computer, Inc., registered in the U.S. and other countries. Sony is a registered trademark of Sony Corporation. Windows is a registered trademark of Microsoft Corporation in the United States and other countries. Other brand and product names are the trademarks or registered trademarks of their respective owners and manufacturers.



USA 520.889.9933
Asia Pacific +65.6841.2094

France +33.1.60.86.03.65
Germany +49.89.660.779.3

Japan +81.3.5639.2731
UK +44.1628.890858

 **PHOTOMETRICS®**
www.photometrics.com
info@photometrics.com