



LVDS PCIe Card

March 23, 2015

Photometrics released its 64-bit only compatible LVDS PCIe card for current Photometrics LVDS-based cameras.

The PCI-Express LVDS card extends the lifetime of existing Photometrics cameras by providing full support for newer and more current computer systems

The 64-bit LVDS card requires availability of a PCI-Express **1X slot** in the computer system.

The PCI-Express LVDS interface card provides 64-bit support for current Photometrics LVDS-based cameras. Additional information is available on the Photometrics website:

- [PCIe LVDS Card Datasheet](#)
- [Photometrics Computer Recommendations](#)
- The PVCAM camera driver is required. [Download a free copy of PVCAM 2.9.3 or higher.](#)

Quote:

Part Number: PCIE-CARD-LVDS

Description: PCI-Express LVDS Card for Photometrics Cameras.

Contact Photometrics Support for Pricing:

[Online Request Form](#)

Phone: +1 800.874.9789

Email: support@photometrics.com

Order:

- [Complete and submit the Online Request Form](#) to place and process an LVDS PCIe card order.

HIGH PERFORMANCE EMCCD & CCD CAMERAS FOR LIFE SCIENCES