



Supported IEEE FireWire and USB Cards

IEEE FireWire provides sustained data rates without requiring a computer host for interconnection between peripherals. These are the IEEE-400 cards that are recommended by Photometrics:

**FireBoard Blue-e™ 1394a OHCI
PCI-Express Adapter**
(Part number: 01-IEEE-PCIe)



**FireBoard Blue™ 1394a
PCI Adapter**
(Part number: 01-IEEE-PCI)



USB provides high speed data transfer. This is the USB card recommended by Photometrics:

2 Port PCI Express (PCIe) SuperSpeed USB 3.0 Card Adapter with UASP - LP4 Power
(Part number: 01-USB3-PCIE)



HIGH PERFORMANCE EMCCD & CCD CAMERAS FOR LIFE SCIENCES

Note:

Photometrics recommends that power be connected from the PC power supply to the card power port via a Y- Cable and or sata to IDE adaptor. Photometrics recommends also ordering this product for use with the card.

(Part numbers: 11-0127 & 16-0150)



If you have any questions, please contact Photometrics Support:

- [Online Request Form](#)
- Phone: +1 800.874.9789
- Email: support@photometrics.com

To Order:

[Complete and submit the Online Request Form](#) to place and process your order.

HIGH PERFORMANCE EMCCD & CCD CAMERAS FOR LIFE SCIENCES