



## Updating the FireWire 1394 Bridge Firmware

1. Close any programs using the camera. If any applications are open, you will be unable to update the firmware.
2. Extract firmware.zip to a location of your choice (eg. Desktop).




3. Run "PMFWUtil" from the Photometrics folder in the Start Menu.







4. Make sure the "Photometrics 1394 Bridge" is selected for the Node.



5. Click "Firmware" under the Load section



6. Select the Firmware file which was unzipped in step 2.



7. The firmware will now update. Once it is completed, exit the program.

